




ANTERIOR
IP CE 600171-A2

Saw Blades for the
PLUS SOLUTION KNEE FAMILY

 **smith&nephew**
PLUS SOLUTION
KNEE FAMILY
Complete Knee System



PLUS SOLUTION
KNEE FAMILY

Saw Blades for the PLUS SOLUTION KNEE FAMILY

The saw blades from Smith & Nephew have been specially developed for the PLUS SOLUTION KNEE FAMILY together with Gebr. Brasseler/Komet Medical. The tothing patented by Gebr. Brasseler guarantees a precise cut and mini-mises the risk of thermal necrosis. Torsion-proof stainless steel allows flat cutting and makes a major contribution towards the fit of our implants.

We are offering a sterile version of the blade for single use. The packaging concept meets the latest sterile-packaging requirements prevailing throughout the world.

Depth of cutting edge: 90 mm

Linvatec Hall

Hall 3/4 – Versipower (Plus) –
Power Pro – MPower 6300

DeSoutter

new – Multi Drive KDX-600

Linvatec 3M

Maxi Driver

Synthes/AO
















Compact Air Drive – Power Drive
Air Power Drive – Battery Power Line

Stryker

Stryker-System 2000/4/5 – EHD

Aesculap

Acculan – Comfort GB432

	Width of cutting edge	Thickness 1.0 mm		Thickness 1.27 mm	
		Art. No. Plus	Art. No. S&N	Art. No. Plus	Art. No. S&N
	13 mm	260084	75006152	260090	75006158
 new	19 mm	–	–	KMS 1913.01S	75023012
	22 mm	260074	75006146	260095	75006163
	13 mm	260079	75006148	260091	75006159
	22 mm	260069	75006142	260093	75006161
	13 mm	260080	75006149	260088	75006156
	22 mm	260070	75006143	260094	75006162
	13 mm	260081	75006150	260086	75006154
 new	19 mm	–	–	KMS 1913.51S	75023011
	22 mm	260071	75006144	260097	75006165
 new Eazy Cut (non-sterile)	22 mm	–	–	EZ 2213.57F	75023015
	13 mm	260083	75006151	260089	75006157
 new	19 mm	–	–	KMS 1913.61S	75023010
	22 mm	260073	75006145	260096	75006164
	13 mm	260085	75006153	260087	75006155
 new	19 mm	–	–	KMS 1913.81S	75023013
	22 mm	260075	75006147	260092	75006160
 new Eazy Cut (non-sterile)	22 mm	–	–	EZ 2213.87F	75023014

Distributor

Smith & Nephew Orthopaedics AG
Erlenstrasse 4a
6343 Rotkreuz
Switzerland

For further information please
contact our local sales office.
www.smith-nephew.com