




**KOMET MEDICAL**  
surgical power accessories



# Große Knochenchirurgie - Sägeblätter für Aesculap®-Antriebe

*Big bone surgery - Saw blades for Aesculap®-power systems*

Blattstärke Material thickness (mm)	<b>KMS</b> sterile (single use) 10 p. box		<b>EZblade®</b> Premium   sterile (reusable) 1 p. box		Schnitttiefe Cutting depth	Schnittbreite Cutting width
	Acculan alt/old	3Ti	Acculan alt/old	3Ti		
1,00	KMS1310.88R	KMS1310.M81	EZ1310.M80	EZ1310.M81	90 mm	13 mm
1,00	KMS2210.86F	KMS2210.M81	EZ2210.M80	EZ2210.M81	90 mm	22 mm
1,00	KMS2510.811	KMS2510.M81	-	-	90 mm	25 mm
1,19	KMS1312.89J	-	-	-	90 mm	13 mm
1,19	-	-	EZ2212.M80	EZ2212.M81	90 mm	22 mm
1,19	KMS2512.85J	KMS2512.M81	-	-	90 mm	25 mm
1,27	KMS1313.807	KMS1313.M81	EZ1313.M80	EZ1313.M81	90 mm	13 mm
1,27	KMS1713.L80	-	-	-	100 mm	17 mm
1,27	KMS1913.81S	KMS1913.M81	EZ1913.M80	EZ1913.M81	90 mm	19 mm
1,27	KMS2213.87F	KMS2213.M81	EZ2213.M80	EZ2213.M81	90 mm	22 mm
1,27	KMS2513.812	KMS2513.M81	-	-	90 mm	25 mm
1,37	-	-	EZ1314.M80	EZ1314.M81	90 mm	13 mm
1,37	-	-	EZ2214.M80	EZ2214.M81	90 mm	22 mm
1,37	KMS2514.85T	KMS2514.M81	-	-	90 mm	25 mm
1,47	-	-	EZ2215.M80	EZ2215.M81	90 mm	22 mm
1,47	KMS2515.813	KMS2515.M81	-	-	90 mm	25 mm



Aesculap ist die Handelsmarke von Aesculap AG & Co. KG. KOMET MEDICAL ist nicht affiliert mit Aesculap.  
Aesculap is a trademark of Aesculap AG & Co. KG. KOMET MEDICAL is not affiliated with Aesculap.



# Kompatible Sägeblätter für Acculan 3Ti

compatible saw blades for Acculan 3Ti

101

**REF**  
Materialdicke  
Material Thickness

<b>KMS1910.S81</b>	<b>1,00 mm</b> 0.039 "				
--------------------	---------------------------	--	--	--	--

101

**REF**  
Materialdicke  
Material Thickness

<b>KMS2510.S81</b>	<b>1,00 mm</b> 0.039 "				
--------------------	---------------------------	--	--	--	--

**REF**  
Materialdicke  
Material Thickness

<b>KM51910.C81</b>	<b>0,40 mm</b> 0,016 "	Schnittstärke Cutting Thickness	<b>0,60 mm</b> 0,024 "		
--------------------	---------------------------	------------------------------------	---------------------------	--	--

**REF**  
Materialdicke  
Material Thickness

<b>KM51915.C81</b>	<b>0,40 mm</b> 0,016 "	Schnittstärke Cutting Thickness	<b>0,60 mm</b> 0,024 "		
--------------------	---------------------------	------------------------------------	---------------------------	--	--

**REF**  
Materialdicke  
Material Thickness

<b>KM51917.S81</b>	<b>0,60 mm</b> 0,024 "	Schnittstärke Cutting Thickness	<b>0,80 mm</b> 0,031 "		
--------------------	---------------------------	------------------------------------	---------------------------	--	--

**REF**  
Materialdicke  
Material Thickness

<b>KM51925.C81</b>	<b>0,40 mm</b> 0,016 "	Schnittstärke Cutting Thickness	<b>0,60 mm</b> 0,024 "		
--------------------	---------------------------	------------------------------------	---------------------------	--	--