




Komet Medical.

















KMS saw blades for joint replacement surgery.



The broadly diversified range of KMS saw blades is suitable for all commonly used power systems.



Explore our Catalogue here.

Widths	Cutting Thicknesses	Connections	Compatible Connections/Hub Styles
8 mm 	0,89 mm	 Aesculap 3Ti	Acculan 3Ti (GA673, GA668), Acculan 4 (GA331, GA336)
10 mm 			
13 mm* 	0,95 mm	 DePuySynthes/ Zimmer Biomet	DePuy Synthes Compact Air Drive II + Power Drive (511.800, 511.801), Air Power Line (511.610), Battery Power Line (530.610), Battery Power Line II (530.710), TRS (05.001.224, 05.001.223, 05.001.240)
17 mm 	1,00 mm		
19 mm 	1,19 mm	 ConMed/ Linvatec/Hall	M Power (PRO6300), M Power 2 (PRO6350M), PowerPro Max (PRO5300M), PowerPro + PowerPro Electric II (PRO6125), PowerPro Pneumatic (PRO6175), Hall Serie 4 (5067-002), Hall Serie 3 (5044-002), Versipower Plus (5071-002), Hall 50 (PRO7300B)
22 mm* 	1,27 mm		
22 mm* 	1,37 mm	 MicroAire	Series 7000
25 mm* 	1,47 mm		
25 mm* 	1,47 mm	 Stryker	System EHD (296-98), System 2000 (2108), System 4 (4108), System 5 (4208), System 6 (6208), System 7 (7208), System 8 (8208)
25 mm* 	1,47 mm		 DeSoutter

*available in parallel and flared