


# Komet Medical<sup>®</sup> twist drills for large bone surgery.

Komet Medical has extended its range of twist drills by adding a new product line. Apart from traditional twist drills with round shank and AO quick coupling, cannulated twist drills are now also part of our range.



SC 28... (cannulated)	(Ø = mm)	(Ø = mm) Cannulation	# FL	HL = mm	Effective length	Connection	REF
	2,0	1,2	3	35	50	AO	SC2805.098
	2,7	1,4	4	50	130	AO	SC2803.098
	3,2	1,7	4	50	140	AO	SC2801.098
	3,5	1,4	4	50	130	AO	SC2804.098
	4,5	1,7	4	50	140	AO	SC2802.098

# Range | Twist drills



S27...	(Ø = mm)	# FL	HL = mm	Effective length	Connection	REF
	1,1	2	16	30	AO	S2709.098
	1,5	2	30	55	AO	S2710.098
	1,5	2	30	80	AO	S2711.098
	1,5	2	30	95	AO	S2712.098
	1,8	2	30	70	AO	S2713.098
	2,0	2	50	70	AO	S2714.098
	2,0	2	50	95	AO	S2715.098
	2,4	2	50	70	AO	S2716.098
	2,5	2	50	150	AO	S2717.098 (golden)
	2,5	2	50	80	AO	S2718.098 (golden)
	2,7	2	50	70	AO	S2719.098
	2,7	2	50	95	AO	S2720.098
	3,2	2	50	165	AO	S2721.098
	3,2	2	50	115	AO	S2722.098
	3,2	2	50	195	AO	S2723.098
	3,5	2	50	80	AO	S2724.098
	3,5	2	50	165	AO	S2725.098
	4,0	2	50	195	AO	S2726.098
	4,5	2	50	115	AO	S2727.098
	4,5	2	50	165	AO	S2728.098
	6,0	2	50	165	AO	S2729.098



S27...	(Ø = mm)	# FL	HL = mm	Effective length	Connection	REF
	1,1	2	20	20	Jakobs	S2730.501
	1,5	2	40	40	Jakobs	S2731.501
	2,0	2	50	55	Jakobs	S2732.501
	2,5	2	50	115	Jakobs	S2733.501
	2,5	2	50	55	Jakobs	S2734.501
	2,5	2	50	75	Jakobs	S2735.501
	2,7	2	50	55	Jakobs	S2736.501
	2,8	2	50	125	Jakobs	S2785.501
	3,2	2	60	100	Jakobs	S2737.502
	3,2	2	60	150	Jakobs	S2738.502
	3,2	2	50	108	Jakobs	S2788.502
	3,5	2	60	70	Jakobs	S2739.502
	3,5	2	50	90	Jakobs	S2724.501
	3,5	2	60	150	Jakobs	S2740.502
	4,5	2	80	100	Jakobs	S2741.502
	4,5	2	80	150	Jakobs	S2742.502